


# The NEW SAPPHIRE 4.0



New LEDs + New Lenses + New Control, contactless control and with wireless remote control


More efficient heat management with the latest technologies and materials

New ergonomic design with smaller housing cross-section and diameter

Parallel suspension and 10kg less weight ensure perfect balance and positioning







### Advantages

- Unparalleled product quality and Know How
- Extreme long life time
- Extended warranty
- Maintenance free
- Extreme low running costs
- Environmental friendly production and materials
- **100% Made in Germany**

### Standard Features

- Stepless dimming
- Focusable light field
- Control of colour temperature
- Unique contactless and sterile control of all functions via sensor technology
- Wireless remote control (optional)
- Several HD- and 4K camera options

Light Intensity	160.000lx at 1m distance
Dimming Range	5% - 100%
Light Field Diameter	200mm - 380mm
Colour Temperature	3000K - 5000K stepless adjustable
Colour Rendering Index	>97
Colour Rendering R9	>97
UV-Radiation	None
IR-Radiation	None
Illumination Depth	1530mm
Supply Voltage	24V / DC
Power Consumption	Max 160W
Mains Power	110 - 230V adjustable
Life Time	>60.000h
Size of Light Head	590 x 590mm

Tolerance  $\pm 10\%$ ; technical specifications are subject to change

