

RUBY-plus



Advantages

- Unbeatable product quality and Know How
- Extreme long life time
- Extended warranty
- Maintenance free
- Extreme low running costs
- Environmental friendly production and materials
- **100% Made in Germany**

Standard Features

- Stepless dimming
- Control of colour temperature
- Unique touch less and sterile control of all functions via sensor technology
- Endo-light
- Sterile handle

Light Intensity	85.000lx at 1m distance
Dimming Range	5% - 100%
Light Field Diameter	Ø200mm at 1m distance
Colour Temperature	3.000K - 5.000K adjustable
Colour Rendering Index	>96
Colour Rendering R9	>96
UV-Radiation	None
IR-Radiation	None
Illumination Depth	800mm
Supply Voltage	24V / DC
Power Consumption	max. 40W
Mains Power	110V - 230V, adjustable
Life Time	>50.000h
Size of Light Head	300 x 300mm

