


ZIRCON



Advantages

- Unparalleled product quality and Know How
- Extreme long life time
- Extended warranty
- Maintenance free
- Extreme low running costs
- Environmental friendly production and materials
- **100% Made in Germany**

Standard Features

- Stepless dimming
- Control of colour temperature
- Unique touch less and sterile control of all functions via sensor technology
- Wireless remote control (optional)
- Several HD- and standard camera options
- Endo-light

Light Intensity	100.000lx - 120.000lx at 1m
Dimming Range	5% - 100%
Light Field Diameter	280mm
Colour Temperature	3000K - 5000K adjustable
Colour Rendering Index	>96
Colour Rendering R9	>96
UV-Radiation	None
IR-Radiation	None
Illumination Depth	1100mm
Supply Voltage	24V / DC
Power Consumption	Max 100W
Mains Power	110 - 230V adjustable
Life Time	>50.000h
Size of Light Head	570 x 570mm

